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SOCOREX
SWISS

STEPPER™

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Sterilization and Purity Certificate – Ecostep *bioproof™*

Socorex® certifies that Ecostep *bioproof™* syringes, intended for use with the Stepper™ repeater pipette are sterilized and free of detectable human DNA, DNase, RNase and Pyrogen (endotoxins).

Quality controls performed on each lot by independent laboratories according to procedures below.

Sterility SAL (sterility assurance level): 10⁻⁶

Sterilization ref: CSIR S.A. SOP ref 7.5.1-13, page 1
Method: Gas mixture 90% E2O, 10% CO₂, initial vacuum 50 mbar, temperature 50°C, 90, humidity 90%. Theoretical calculated gas concentration 77.2% E2O / L, exposure time 4h, 5 rinsing steps
Sterility test ref: CSIR S.A. SOP 7.5.1-30, ISO 11138-1 and -2
Method: Minimum 18 self-contained biological indicators Bacillus Alpbreus (Subtilis var. niger), ATCC No. 35762 at 10⁶ concentration, 2 days incubation at 37°C
Validation: No growth detected

Human Deoxyribonucleic Acid (DNA) < 2 pg

Test ref: Socorex Research S.A., SOP ref. SAM5003 v1, 2014-01
Method: Amplification by PCR of "Alu" genomic area in human DNA. Migration of PCR product on agarose gel. Test performed on Ecostep™ syringes rinsed with DNA-free water.
 1) Syringe from tested lot 2) Syringe added with 1 µg of Human DNA.
 3) Negative control 4) Positive control
Validation: No DNA amplification in samples 1) and 3), detected amplification in samples 2) and 4)

Deoxyribonuclease (DNase) < 10⁷ Kunitz Units

Test ref: Socorex Research S.A., SOP ref. SAM5001 v1, 2014-01
Method: Incubation on agarose gel of DNA molecular scale. Test performed on Ecostep™ syringes.
 1) Syringe from tested lot 2) Syringe added with a 10⁷ K Unit of DNase
 3) Negative control 4) Positive control
Validation: No degradation on agarose gel of the DNA molecular scale in samples 1) and 3), degradation in samples 2) and 4)

Ribonuclease (RNase) < 10⁶ Kunitz Units

Test ref: Socorex Research S.A., SOP ref. SAM5002 v1, 2014-01
Method: Incubation on agarose gel of RNA molecular scale. Test performed on Ecostep™ syringes.
 1) Syringe from tested lot 2) Syringe added with a 10⁶ K Unit of RNase
 3) Negative control 4) Positive control
Validation: No degradation on agarose gel of RNA molecular scale in samples 1) and 3), degradation in samples 2) and 4)

Pyrogen (endotoxins) < 0.05 IU or EU/mL, < 0.5 EU or EU item tested

Test ref: LAL chromogenic method, European Pharmacopoeia 8th edition (2014), chapter 2.6.14, and United States Pharmacopoeia 37 NF 32 (2014), chapter 85
Method: Preparation of a standard curve from 5 IU (or EU/mL) to 5 · 10⁻⁶ IU (or EU/mL). Bacterial endotoxin rates determined using spectrophotometric measure at 600 nm.
 Test performed on Ecostep™ syringes rinsed with Pyrogen-free water.
 1) Syringe from tested lot 2) Syringe added with 0.5 IU (or EU/mL) of endotoxin
 3) Negative control
Validation: No selection in samples 1) and 3), detection in sample 2)

LOT SPECIFIC CERTIFICATE AVAILABLE FROM: SOCOREX@SOCOREX.COM

SOCOREX
 socorex@socorex.com www.socorex.com Ecolab, September 2020
 Jean-Marc Ammann
 Quality System &
 Regulatory Affairs Manager

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Copyright © 2020 by Socorex Isba SA, Switzerland.