CALIBRA[®] 822

DIGITAL MICROPIPETTE


DIGITAL MIKROPIPETTE

MICROPIPETTE DIGITALE

OPERATING INSTRUCTIONS

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MODE D'EMPLOI

SAFETY INSTRUCTIONS

- Read operating instructions carefully; keep available for future reference
- Observe manufacturer's instructions and guidelines for operation and maintenance
- Refer to and follow regulations about handling of potentially hazardous reagents
- Before use, check instrument for good working conditions
- Do not use harsh chemicals (such as acetone) to clean instrument
- This product should be used only for its intended purpose
- Mind possible hand-fatigue during serial pipetting and its medical consequences (such as repetitive strain injuries RSI)

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DESCRIPTION

The Calibra[®] 822 is an air displacement adjustable volume micropipette intended for the accurate, precise and easy pipetting of volumes from 10 to 1000 µL. Its unique, patented Solid Calibration[®] Twin Cam system provides for both digital entry and digital display of the volume. Materials have been carefully selected to ensure durability and optimal resistance to aggressive media.

The unit is autoclavable at 121°C fully assembled.

Other major features include instant volume setting, integral tip ejector activated by separate push button, slim shaft for narrow vials and attractive colour coding. Calibration stability owed to precalibrated Twin Cam System. However in-lab calibration easily performed when required.

VOLUME SETTING (Fig. 1)

Two cylindrical cams fitted with pre-calibrated steps rotate on their axis. The large cam shows the figures on the left hand side of the window, the small one on its right hand side.

The volume setting as described below and illustrated in Fig. 1 is fast, easy and glove friendly.

1. Push plunger button with thumb and pointer to the first stop and lock it by a clockwise quarter of a turn.
- 2a. Turn setting wheel counter-clockwise until the desired volume appear at the left of the window.
- 2b. If required, pull out setting wheel (position clicks-in) and repeat point 2a for the fine tuning.
3. Push setting wheel back to the initial position and unlock plunger button.

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Volumen einstellen

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HOLDING MICROPIPETTE (Fig. 2)

The ergonomic shape of the Calibra® greatly reduces operator fatigue. Place fingerrest on the phalanx of forefinger with ejector knob turned towards yourself. Ease of access of both push button and ejector knob allows natural and comfortable thumb movement.

TIPS AND FILTERS

Disposable Qualitix® pipette tips with or without filter were carefully designed to provide the best pipette and tip combination for all types of applications. See p. 4 for ordering information. Detailed flyer on request.

Rinsing every new tip once with the pipetting solution (fill – empty – expulse) will improve precision.

Regular tips are autoclavable at 121°C. Filter tips are either autoclavable or supplied in a sterile version.

TIGHTNESS

To work correctly, the air chamber of the pipette has to be tight. Leaks can be detected by observing the tip: a droplet falling from the extremity signals a leak in the system.

Following spots have to be watched for leakage:

1. Matching of the tip on nozzle
2. O-ring providing tightness to the plunger
3. Damaged nozzle

Leakage will usually be eliminated by changing the tip. If the problem remains, change the O-ring. The wear of the O-ring is very slow under normal use, somewhat faster under frequent autoclaving. Sterilization therefore requires higher frequency tightness control. A damaged nozzle should be changed.


Fig. 2

PIPETTING OPERATION (Fig. 3)

Forward pipetting mode:

- Phase 1** Depress plunger button to the intermediate stop ending calibrated stroke before immersing tip into appropriate reagent container.
- Phase 2** Immerse tip vertically 2-3 mm into liquid and smoothly release plunger button. Withdraw tip from liquid vertically after 2-3 seconds.
- Phase 3** Place tip against wall of receiving vessel, depress plunger button to first stop and let liquid out. Depress button to second stop: the excess air volume produced by depressing the plunger beyond the calibrated stroke blows out possible fluid particles which may cling to the tip.

Tip ejection Eject used tip by pressing ejector (E) as shown in Fig. 3. Change tip per recommendations in previous section.

Reverse pipetting mode also possible with the Calibra® pipettes.

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CALIBRATION

Each Calibra® has been calibrated permanently at the factory and individually controlled. The Q.C. certificate accompanying the instrument includes control data, serial number and operator identification. Calibration is carried out at a constant temperature (+/- 0.5°C) comprised between 20 and 25°C with Aqua bidest. according to ISO 8655.

Under normal conditions, a Solid Calibration pipette does not need to be recalibrated. After replacing a part, which could influence the measure, the calibration should be checked, preferably by gravimetry according to the above mentioned norm; an analytical balance is needed. Ten measurements should be performed after checking the seal tightness of the instrument.

In-lab calibration (Fig. 4): Whenever required, for instance when pipetting viscous or non-aqueous liquids, calibration adjustments can be performed by laboratory personnel. Reach calibration screw c) located inside plunger button by removing comfort pad a). By means of allen key b), delicately turn calibration screw clockwise to decrease volume, and vice versa, using values indicated below. Check resulting volume after each fine tuning. Cover push-button with comfort pad again.

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
ECHELLE DE CALIBRATION

MAINTENANCE

The Calibra® 822 was designed to sustain long-term operation with minimal maintenance and service requirements. Observing the following points will ensure best instrument performance at all times:

- **Clean and dry air channel and plunger** if liquid accidentally entered into the instrument. To do so, unscrew lower part of the pipette, clean air channel with rod supplied if necessary, use alcohol or appropriate detergent, rinse air channel and plunger, allow parts to dry before greasing and reassembling.
- **Replace O-ring** or nozzle as may be required by lack of tightness (refer to tightness section above).

Clean, grease and check instrument once a year or more depending upon frequency of use.

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STERILIZATION

The unit has been designed for repeated sterilisation in the autoclave at 121°C (20 minutes), fully assembled. Before autoclaving, set pipette to the largest possible volume. Place instrument horizontally in the autoclave avoiding any direct contact with metal.

Allow instrument to dry and cool down before use. Tighten lower assembly nut if loose. If plunger movement is not smooth, open the pipette and slightly lubricate plunger with silicone grease supplied with the instrument.

Check tightness and accuracy of the pipette regularly but at least after 50 autoclaving cycles.

A slight change of material colours may appear after long term repeated autoclaving. Metal nozzle of 2 and 10 µL models can be briefly flamed for microbiological applications. Correct autoclaving and resulting sterility are the responsibility of the user.

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WARRANTY

Your Calibra® digital pipette is guaranteed against any material or manufacturing defects for the period of time specified in its QC certificate. Damages due to non-respect of manufacturer's instructions, safety precautions or autoclaving conditions, as well as material colour alteration are excluded from the warranty. Repair and replacement of parts do not extend warranty time. Should regular maintenance not eliminate a detected defect, return the instrument to the dealer from whom it was purchased after obtaining return authorization. **Decontaminate instrument prior to returning it.**

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PERFORMANCE DATA AND ORDERING INFORMATION

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* Values obtained with forward pipetting mode, at constant temperature (+/-0.5°C) comprised between 20 and 25°C with Aqua bidest according to ISO 8655. Viscous or volatile liquids may generate deviations greater than shown in above figures.

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