# CALIBRA<sup>®</sup> 832

## DIGITAL MACROPIPETTE


## DIGITAL MAKROPIPETTE

## MACROPIPETTE DIGITALE

### OPERATING INSTRUCTIONS

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### MODE D'EMPLOI

#### SAFETY INSTRUCTIONS

- Read operating instructions carefully; keep available for future reference
- Observe manufacturer's instructions and guidelines for operation and maintenance
- Refer to and follow regulations about handling of potentially hazardous reagents
- Before use, check instrument for good working conditions
- Do not use harsh chemicals (such as acetone) to clean instrument
- This product should be used only for its intended purpose
- Mind possible hand-fatigue during serial pipetting and its medical consequences (such as repetitive strain injuries RSI).

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#### DESCRIPTION

The Calibra<sup>®</sup> 832 is an air displacement, adjustable volume macropipette intended for the accurate, precise and easy pipetting in the 0.2 – 2 mL and 1 – 10 mL ranges. Its unique, patented Solid Calibration<sup>®</sup> Twin Cam system provides for both digital entry and digital display of the volume.

The materials of the Calibra<sup>®</sup> 832 have been carefully selected to resist chemical corrosion due to fumes and emanations of aggressive media. Cellulose nozzle filters are available for added safety. The unit is autoclavable at 121°C fully assembled.

The liquid sample is drawn into a disposable polypropylene tip (also Pasteur pipette on the 0.2 – 2 mL model, adapter supplied with the instrument) without coming in contact with the instrument. Fully integrated in the handpiece, the tip ejector releases the tip automatically. Each instrument is individually calibrated and comes along with its own control certificate.

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#### VOLUME SETTING (Fig. 1)

Two cylindrical cams fitted with pre-calibrated steps rotate on their axis. The large cam shows the figures on the left hand side of the window, the small one on its right hand side.

The volume setting as described below and illustrated in Fig. 1 is fast, easy and glove friendly.

1. Push plunger button with thumb and pointer to the first stop and lock it by a clockwise quarter of a turn.
- 2a. Turn setting-wheel counter-clockwise until the expected key volume appears at the left of the window.
- 2b. If required, pull out setting wheel (position clicks-in) and repeat point 2a for the fine tuning.
3. Push setting wheel back to the initial position and unlock plunger button.

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Abb. 1. Volumen einstellen

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## HOLDING THE MACROPIPETTE

The ergonomic shape of the Calibra® greatly reduces operator fatigue. Place fingerrest on the phalanx of forefinger with ejector knob turned towards yourself. Ease of access of both push-button and ejector knob allows natural and comfortable thumb movement.


Fig. 2

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
Fig. 4

## TIPS AND FILTERS

For consistent results and proper operation of the ejector, use original 2 mL, respectively 10 mL macro tips. Made of virgin, dye and metal free polypropylene (PP), they are autoclavable at 121°C. A special adapter supplied with the 0.2 – 2 mL model accommodates standard Pasteur pipettes (ext. Ø 6.5 – 7.2 mm) as well. Glass Pasteur pipettes are of advantage in handling PP affecting solvents. Refer to page 4 for details.

Apply new tip or Pasteur pipette as often as deemed necessary in regard with changes of sample nature, droplet accumulation after blow out, or foam formation. Enhance precision by filling new tip or Pasteur pipette once before sampling.

The use of optional nozzle filters drastically reduces the risks of barrel and plunger contamination. Their use is particularly recommended with large volumes and/or when handling biohazardous, radioactive or highly aggressive liquids. Change filter if wet or contaminated. Filters must not be autoclaved.

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## PIPETTING OPERATION (Fig. 2)

Forward pipetting mode:

**Phase 1** Depress plunger button to the intermediate stop ending calibrated stroke before immersing tip into appropriate reagent container.

**Phase 2** Immerse tip vertically 5 – 10 mm into liquid and smoothly release plunger button. Withdraw tip from liquid vertically after 2-3 seconds, sliding it along container wall.

**Phase 3** Place tip against wall of receiving vessel, depress plunger button to first stop and let liquid out. Depress button to second stop: the excess air volume produced by depressing the plunger beyond the calibrated stroke blows out possible fluid particles which may cling to the tip.

**Tip ejection** Eject used tip by pressing ejector (E)

as shown in Fig. 2. Change tip per recommendations in previous section.

Reverse pipetting mode also possible with the Calibra® macro pipette.

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## TIGHTNESS

To work correctly, the air chamber of the pipette has to be tight. Leaks can be detected by observing the tip: a droplet forming and falling from the extremity signals a leak in the system.

Following spots have to be watched for leakage:

1. Matching of the tip on nozzle.
2. X-ring providing tightness on the plunger.
3. Pasteur pipette adaptor's O-ring providing tightness.
4. Damaged barrel inner wall or nozzle.

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Leakage will be usually eliminated by changing the tip. If the problem remains, change the X-ring. The wear of the X-ring or O-ring is very slow under normal use, somewhat faster under frequent autoclaving. Therefore, autoclaving requires higher frequency tightness control. A damaged barrel or nozzle should be changed.

## CALIBRATION

Each Calibra® has been calibrated permanently at the factory and individually controlled. The Q.C. certificate accompanying the instrument includes control data, serial number and operator identification. Calibration is carried out with the protection nozzle filter, at a constant temperature (+/- 0,5°C) comprised between 20 and 25°C with Aqua bidest. according to EN ISO 8655.

Under normal conditions, a Solid Calibration® pipette does not need to be recalibrated. After replacing a part, which could influence the measure, the calibration should be checked, preferably by gravimetry according to the above mentioned norm; an analytical balance is needed. Ten measurements should be performed after checking the seal tightness of the instrument.

**In-lab calibration (Fig. 3):** Whenever required, for instance when pipetting viscous or non-aqueous liquids, calibration adjustments can be performed by laboratory personnel. Reach calibration screw c) located inside plunger button by removing comfort pad a). By means of allen key b), delicately turn calibration screw clockwise to decrease volume, and vice versa using data shown in Fig. 4. Check resulting volumes after each fine tuning. Cover push-button with comfort pad again.

## MAINTENANCE

The Calibra® 832 pipettes are in conformity with Council directive IVD 98/79 EEC (in vitro diagnostic medical devices). Consequently, any change on the pipettes (i.e. parts, calibration) is done under responsibility of the user. The models were designed to sustain long-term operation with minimal maintenance and service requirements. Observing the following points will ensure best instrument performance at all times:

- **Clean and dry barrel and plunger** if liquid accidentally entered into the instrument. Dry sample residues may damage plunger seal ring and barrel wall. Rinsing barrel and plunger after handling highly aggressive fluids will lengthen instrument life (N.B.: the use of protection nozzle filters is recommended with such liquids and/or when working with large volumes). Do not use any tool nor a rough brush when cleaning inner barrel wall. Scratches would reduce air seal tightness and require barrel replacement.
- **Lubricate barrel wall** with supplied silicon grease from time to time; a detergent cleaned or new barrel is to be lubricated prior to installation.
- **Replace plunger ring (X-ring)** if damaged. Remove ring by squeezing it side wise between thumb and pointer, pull it out with nails or paper clip. Install new ring.
- **Replace Pasteur pipette adapter** if pipette fits loose or worn out X-rings cause a lack of tightness.

Clean, grease and check instrument once a year or more depending upon frequency of use.

## STERILIZATION

The unit has been designed for repeated sterilization in the autoclave at 121°C (20 minutes) fully assembled. Before autoclaving, set pipette to the largest possible volume and remove nozzle filter. Place instrument horizontally in the autoclave avoiding any direct contact with metal.

Allow instrument to dry and cool down before use. Tighten cylinder screw of 10 mL model if loose. If plunger movement is not smooth, open the pipette and slightly lubricate plunger and barrel with silicone grease supplied with the instrument. Change nozzle protection filter.

Check tightness and performance after first autoclaving cycle then regularly but at least after 50 autoclaving cycles.

A slight change of material colours may appear after long term repeated autoclaving.

Correct autoclaving and resulting sterility are the responsibility of the user.

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## WARRANTY

Your Calibra® digital macropipette is guaranteed against any material or manufacturing defects for the period of time specified in its QC certificate. Damages due to non-respect of manufacturer's instructions, safety precautions or autoclaving conditions, as well as material colour alteration are excluded from the warranty. Repair and replacement of parts do not extend warranty time. Should regular maintenance not eliminate a detected defect, return the instrument to the dealer from whom it was purchased after obtaining return authorisation. **Decontaminate instrument prior to returning it.**

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## PERFORMANCE DATA AND ORDERING INFORMATION

Volume	Division	Inaccuracy (E%)			Imprecision (CV%)			Cat. No.
		Min. vol.	Mid. vol.	Max. vol.	Min. vol.	Mid. vol.	Max. vol.	
0.2 – 2 mL	0.01 mL	< ± 1.5%	< ± 1.0%	< ± 0.5%	< 0.5%	< 0.3%	< 0.2%	832.02
1 – 10 mL	0.1 mL	< ± 1.5%	< ± 1.0%	< ± 0.5%	< 0.3%	< 0.3%	< 0.15%	832.10

\* Values obtained with forward pipetting mode, at constant temperature (+/-0.5°C) comprised between 20 and 25°C with Aqua bidest according to ISO 8655 and nozzle protection filter. Viscous or volatile liquids may generate deviations greater than shown in above figures.

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

**DISASSEMBLING 832.02 (Fig. 5.1 - 5.4)**

1. See Fig. 5.1 Twist out ejector a).
2. See Fig. 5.2 Pull locking catch b) with key c). Catch is not removable.
3. See Fig. 5.3 Disengage assembly d) by pulling hood e) firmly.
4. See Fig. 5.4 Remove universal joint f) and plunger g) if necessary.

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**ASSEMBLING 832.02 (Fig. 5.4 - 5.1)**


1. See Fig. 5.4 Assemble universal joint f) and plunger g).
2. See Fig. 5.3 Engage assembly d) while maintaining hood e) down.
3. See Fig. 5.2 Push in locking catch b) fully.
4. See Fig. 5.1 Click in ejector a).

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**DISASSEMBLING 832.10 (Fig. 6.1 - 6.4)**

1. See Fig. 6.1 Twist out ejector a), unscrew and remove cylinder b).
  2. See Fig. 6.2 Remove universal joint c) and plunger d) if required.
- Note: operations 3 and 4 usually unnecessary
3. See Fig. 6.3 Grasp locking catch e) with key f) and pull it out completely.
  4. See Fig. 6.4 Disengage assembly g) by pulling hood h) firmly.

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**ASSEMBLING 832.10 (Fig. 6.4 - 6.1)**


1. See Fig. 6.4 Engage assembly g) while maintaining hood h) down.
2. See Fig. 6.3 Insert locking catch e) fully; use key f) to let it rotate to clicked in position.
3. See Fig. 6.2 Assemble universal joint c) and plunger d).
4. See Fig. 6.1 Screw in barrel b) and click in ejector a).

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We declare under our sole responsibility that our devices

**CALIBRA®digital pipettes**

are in conformity with the following Council directive

---

**IVD 98/79/EEC**  
(in vitro diagnostic medical devices)

---

Sylvain Christen  
Director, CEO

Jean-Pierre Uldry  
Director, Engineering

Ecublens  
September 3<sup>rd</sup>, 2008

**▲ SOCOREX**

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**SOCOREX ISBA S.A.**