

The two universal micropipette stands fit all brands, ensuring their proper storage at any time. Anti-slip pads guarantee bench stability. Easy to clean.



336

Twister™ universal 336

- 6 positions
- 360° soft rotation
- Instruments easy to pick-up
- Concept fits pipettes of all brands
- Interchangeable discs
- 7 translucent colours



Universal storage solution

No more micropipettes laying on the bench thanks to universal stands. They accommodate Socorex and other brands micropipettes. Simple disassembling and cleaning.

Ordering information


Colour	Packaging	Cat. No.
Twister™ universal 336		
Rubin red	1/ pk	320.336R
Topaz orange	1/ pk	320.336O
Citrine yellow	1/ pk	320.336Y
Emerald green	1/ pk	320.336G
Sapphire blue	1/ pk	320.336B
Quartz grey	1/ pk	320.336Q
Diamond white	1/ pk	320.336W



337

In-line™ universal 337

- 7 positions
- Holds single channel pipettes of all brands
- Space-saving pipette stand
- Sturdy polyamide material
- 5 attractive pastel colours



In-line™ universality

One fits all; the stand holds new and old generations of Socorex and other brands micropipettes.

Ordering information

Colour	Packaging	Cat. No.
In-line™ universal 337		
Light grey	1/ pk	320.337G
Pastel rose	1/ pk	320.337R
Ice blue	1/ pk	320.337B
Mint green	1/ pk	320.337M
Vanilla yellow	1/ pk	320.337Y


Proper micropipette storage prevents damages and extends their life span. Socorex offers a choice of stands adapted to any laboratory configuration.



332

Shelf pipette holder 332

- 2 positions
- Fitting micro, macro and multichannel pipettes



Holder for two instruments

The attractive design accommodates Acura® instruments, including microdispenser pipettes Acura® 865. A self-adhesive strip supplied ensures lasting attachment on any clean surface. Made of polycarbonate material.

Ordering information

Colour	Packaging	Cat. No.
Transparent blue	1 / pk	320.332
Transparent blue	4 / pk	320.332.4

Shelf pipette holder 332

Transparent blue	1 / pk	320.332
Transparent blue	4 / pk	320.332.4



340

Static pipette stand 340

- 3 positions
- Fitting micro, macro and/or multichannel Socorex pipettes
- Easy reach of each instrument



Extended compatibility

Pipette stand 340 also accommodates dilution pipette Acura® 810, positive displacement pipettes Acura® 846 and microdispenser pipettes Acura® 865. Heavy-bottom plate and anti-slip pads guarantee high stability. Made of polyamide material.

Ordering information

Colour	Packaging	Cat. No.
Light grey	1 / pk	320.340

Static pipette stand 340

Light grey	1 / pk	320.340
------------	--------	---------