


# Dosys all-glass syringes

## Model 155

### GENERAL DIMENSIONS

Cat. No.	Volume	Overall length	Cylinder length	Outside diameter
155.0501	1 mL	86	72	9.5
155.0502	2 mL	88	72	11.7
155.0505	5 mL	109	90	14.7
155.0510	10 mL	133	110	18.0
155.0520	20 mL	148	126	23.7
155.0530	30 mL	165	137	27.0
155.0550	50 mL	174	147	32.5
155.05100	100 mL	220	183	39.0
155.05150	150 mL	228	186	48.0
155.05250	250 mL	303	255	48.0

- Borosilicate glass
- Metal Luer Lock Nozzle
- Dimensions are in millimeters



Glass syringes 155\_0912