

manual pipette fillers


Profillette™ 406

Low cost, easy to use pipette filler suitable for graduated glass and plastic pipettes. Rotate wheel to aspirate. Fast release by pressing air inlet valve button, or gradual release by gently turning wheel. One-year warranty.

Assorted package includes one each size.

Product features

- One-hand operation
- Easy filling and dispensing
- Air inlet valve for fast distribution
- Colour coding per size
- Maintenance free



< 2 mL
< 10 mL
< 25 mL

Ordering information - Instruments

Description	Colours	Packaging	Cat. No.
Profillette™ 406, 2 mL	blue	1 / pk	406.002
Profillette™ 406, 10 mL	green	1 / pk	406.010
Profillette™ 406, 25 mL	red	1 / pk	406.025
Set with 2 / 10 / 25 mL	assorted	3 / pk	406.300

Ordering information - Accessories

Description	Packaging	Cat. No.
Spare silicone pipette holder	1 / pk	1.406.01
Replacement nose screw	1 / pk	1.406.02



For more information or for an overview of the complete Socorex product range, visit

www.socorex.com