

# Qualitix® Sterilization paper trays

consumables


Optimal for a large scope of applications in laboratory, medical, dental, clinical, surgery room or veterinary field, the Qualitix® single use paper trays are fully autoclavable by either ETO, Steam or Gamma Ray. Made of high quality 320g/m<sup>2</sup> medical PE coated paper, the anti-skid layer prevents contents from sliding. Ideal for stand-alone use, or inside a Qualitix® sterilization pouch.



## The Qualitix® sterilization paper trays

### Product advantages

- Rigid medical PE coated paper material
- Advantageously replacing stainless steel trays
- One size fits most devices
- Suitable for ETO, steam and gamma sterilization
- Fit for use as such, or in combination with Qualitix® sterilization pouches
- Conformity with ISO 14971, ISO 11135, ISO 17665-06, CE standard



### Convenient use ①

Trays can be used alone or placed inside sterilization pouches, as to preserve sterility after steam, ETO or gamma processing.

### Flexibility ②

A large variety of items can be sterilized, simply placed on rigid paper trays.

### Dimensions ③

Tray size perfectly fits the 190x360 mm dimensions of the sterilization pouch Cat. No. 3210.1936.

Ask for free samples

### Ordering information

Description	Dimensions		Material	Packaging	Cat. No.
Sterilization paper tray	150x265x15 mm	6.0"x10.5"x0.6"	Medical PE coated paper	50 / pk	3220.1526

### Dimensions

A (mm)	B (mm)	C (mm)	D (mm)
150	265	135	250